





Knee Innovations




*GraftLink[®] ACL Reconstruction
Using TightRope[®] II Implant*



*MCL Repair Using the
InternalBrace[™] Technique*





*Meniscal Knee
Scorpion[™] Repair*



*Meniscal Root Repair With PEEK
SwiveLock[®] Anchor*

Knee Scorpion™

Simplified, Single-Portal
Suture Passing and Retrieving



- **Low Profile**
Pass suture in tight spaces with the Knee Scorpion Needle (AR-12990N)
- **Versatility**
Accommodates 2-0 and 0 FiberWire®
- **Efficient**
Pass and retrieve suture in one step through one portal
- **Anatomic**
Angled jaw engineered for the knee joint



Knee Scorpion
AR-12990